

# posterior lateral strut afo


by CUSTOM COMPOSITE MFG, INC



The Posterior Lateral Strut AFO has a similar design to that of the Lateral Strut AFO, but wraps posteriorly, using the calf section of the leg to pick up the toe. It also provides an anteriorly directed force at heel strike to help correct mild recurvatum.

The Posterior Lateral Strut AFO is custom made from a cast of the patient and comes with a foam liner, Velcro® strap and custom foam footbed.

For Casting Instructions, visit [www.cc-mfg.com](http://www.cc-mfg.com) or call us toll free at 866-273-2230.



## indications

Drop foot, MS patients needing mild recurvatum control

## recommended activity level

Low to moderate activity levels

## suggested l-codes

**L-1940** A semi-rigid molded plastic orthosis to hold the foot in neutral position (dorsi-plantar flexion), controls foot position, custom molded from a model of the patient, custom fabricated, includes casting and cast preparation

**L-2755** Addition to lower-extremity orthosis, high strength, lightweight material, all hybrid lamination/prepreg composite, per segment, for custom-fabricated orthosis only

**L-2820** Addition to lower extremity orthosis, soft interface for molded plastic, below knee section